

Citation: Kaku N, Nishimura F, Shigeishi Y, Tachiki R, Sakai H, Sasaki D, et al. (2021) Performance of anti-SARS-CoV-2 antibody testing in asymptomatic or mild COVID-19 patients: A retrospective study in outbreak on a cruise ship. PLoS ONE 16(9): e0257452. https://doi.org/10.1371/journal.pone.0257452

Editor: Etsuro Ito, Waseda University: Waseda Daigaku, JAPAN

Received: June 14, 2021

Accepted: September 1, 2021

Published: September 28, 2021

Copyright: © 2021 Kaku et al. This is an open access article distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited.

Data Availability Statement: Raw data were generated at Nagasaki University Hospital. Derived data supporting the findings of this study are within the paper and its <u>S1</u> and <u>S2</u> tables.

Funding: This work was supported by the Japan Agency for Medical Research and Development (AMED) [grant number JP20he0622041]. The funders had no role in study design, data collection and analysis, decision to publish, or preparation of the manuscript. Cellspect Co., Ltd provided RESEARCH ARTICLE

Performance of anti-SARS-CoV-2 antibody testing in asymptomatic or mild COVID-19 patients: A retrospective study in outbreak on a cruise ship

Norihito Kaku^{1,2¤}*, Fumitaka Nishimura², Yui Shigeishi², Rina Tachiki², Hironori Sakai³, Daisuke Sasaki², Kenji Ota^{1,2}, Kei Sakamoto^{1,2}, Kosuke Kosai^{1,2}, Hiroo Hasegawa^{1,2}, Koichi Izumikawa⁴, Koya Ariyoshi⁵, Hiroshi Mukae⁶, Jiro Yasuda⁷, Kouichi Morita⁸, Shigeru Kohno⁹, Katsunori Yanagihara^{1,2}

 Department of Laboratory Medicine, Nagasaki University Graduate School of Biomedical Sciences, Nagasaki City, Nagasaki, Japan, 2 Department of Laboratory Medicine, Nagasaki University Hospital, Nagasaki City, Nagasaki, Japan, 3 R&D, Cellspect Co., Ltd., Morioka City, Iwate, Japan, 4 Department of Infectious Diseases, Nagasaki University Graduate School of Biomedical Sciences, Nagasaki, Japan,
Department of Clinical Medicine, Institute of Tropical Medicine, Nagasaki University, Nagasaki City, Nagasaki, Japan, 6 Department of Respiratory Medicine, Nagasaki University Graduate School of Biomedical Sciences, Nagasaki City, Nagasaki, Japan, 7 Department of Emerging Infectious Diseases, Institute of Tropical Medicine, Nagasaki University, Nagasaki, Japan, 8 Department of Virology, Institute of Tropical Medicine, Nagasaki University, Nagasaki, Japan, 9 Nagasaki University, Nagasaki City, Nagasaki, Japan

 Current address: Division of Pulmonary and Critical Care Medicine, Department of Internal Medicine, University of Michigan, Ann Arbor, Michigan, United States of America
* norihitk@gmail.com

Abstract

Objectives

A few studies on antibody testing have focused on asymptomatic or mild coronavirus disease 2019 (COVID-19) patients with low initial anti-severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2) antibody responses. Anti-SARS-CoV-2 antibody-testing performance was evaluated using blood samples from asymptomatic or mild COVID-19 patients.

Methods

Blood samples were collected from 143 COVID-19 patients during an outbreak on a cruise ship 3 weeks after diagnosis. Simultaneously, a follow-up SARS-CoV-2 genetic test was performed. Samples stored before the COVID-19 pandemic were also used to evaluate the lateral flow immunochromatographic assay (LFA) and electrochemiluminescence immuno-assay (ECLIA). Titers of anti-SARS-CoV-2 IgM and IgG antibodies against the nucleocapsid and spike proteins were measured using the enzyme-linked immunosorbent assay to confirm which antibodies were influenced on LFA- and ECLIA- false-negative result in crew-member samples.

equipment and reagents for ELISA and support in the form of salaries for author [HS], but did not have any additional role in the study design, data collection and analysis, decision to publish, or preparation of the manuscript. The specific roles of these authors are articulated in the 'author contributions' section.

Competing interests: I have read the journal's policy and the authors of this manuscript have the following competing interests: We used marketed products of Cellspect Co., Ltd for ELISA in this study; Hironori Sakai is an employee of Cellspect Co., Ltd., and performed ELISA in this study; Norihito Kaku and Kosuke Kosai received research funding from Roche Diagnostics K. K.; Katsunori Yanagihara received research funding from Cellspect Co., Ltd., Roche Diagnostics K. K. and, Abbott Japan LLC. These do not alter our adherence to PLOS ONE policies on sharing data and materials.

Results

Sensitivity, specificity, positive-predictive, and negative-predictive values of LFA-detected IgM antibodies were 0.231, 1.000, 1.000, and 0.613, respectively; those of LFA-detected IgG antibodies were 0.483, 0.989, 0.972, and 0.601, respectively; and those of ECLIA-detected total antibodies were 0.783, 1.000, 1.000, and 0.848, respectively. All antibody titers measured using ELISA were significantly lower in blood samples with negative results than in those with positive results in both LFA and ECLIA. In the patients with negative results from the follow-up genetic testing, IgM-, IgG-, and total-antibody positivity rates were 22.9%, 47.6%, and 72.4%, respectively.

Conclusions

These findings suggest that anti-SARS-CoV-2 antibody testing has lower performance in asymptomatic or mild COVID-19 patients than required in the guidelines.

Introduction

The clinical indications for anti-severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2) antibody testing are limited, although many methods have been developed. The detection of anti-SARS-CoV-2 IgG or total antibodies at 3–4 weeks after symptom onset may be useful in determining past infection in selected clinical situations; however, data in this context are limited [1]. On the other hand, anti-SARS-CoV-2 antibody testing may play an essential role in the public health response to coronavirus disease 2019 (COVID-19) and in understanding the outbreak dynamics of the COVID-19 pandemic [2, 3]. However, a previous study reported that initial anti-SARS-CoV-2 antibody responses were lower in asymptomatic or mild COVID-19 patients than in severe COVID-19 patients [4]. Since the percentage of asymptomatic COVID-19 patients in large-sample studies was 43.0%–76.5% [5], it is necessary to elucidate the performance of anti-SARS-CoV-2 antibody testing in asymptomatic or mild COVID-19 patients.

We preserved blood samples of asymptomatic or mild COVID-19 patients collected from a large cruise ship that had experienced a COVID-19 outbreak. The cruise ship anchored at Nagasaki Port from the end of January 2020, and one of the 623 crew members complained of pneumonia symptoms and was diagnosed with COVID-19 in mid-April 2020. Thereafter, SARS-CoV-2 genetic testing was performed on all cruise ship crew members, and SARS-CoV-2 was detected in 144 individuals [6]. However, most patients received care on the ship because they were asymptomatic or did not need inpatient treatment, such as oxygen administration.

In this study, we evaluated the performance of anti-SARS-CoV-2 antibody testing using the lateral flow immunochromatographic assay (LFA) and electrochemiluminescence immunoassay (ECLIA) in blood samples from crew members and those collected from patients between November 2014 and August 2019. In addition, we measured the titers of several anti-SARS--CoV-2 antibodies and compared them between false-negative- and positive-result samples in anti-SARS-CoV-2 antibody testing using the enzyme-linked immunosorbent assay (ELISA) to confirm which antibodies were influenced on LFA- and ECLIA-false-negative results in crewmember samples.

Material and methods

Study design

In this study, asymptomatic or mild COVID-19 was defined as a cruise ship crew member who was not hospitalized when blood samples were collected. Blood samples were collected from 178 crew members in mid-May 2020, 3 weeks after the first SARS-CoV-2 genetic testing. Thirty-four blood samples were excluded from the final analysis due to negative results in the first SARS-CoV-2 genetic testing, and one blood sample was excluded because the patient was hospitalized at the time of sample collection (Fig 1). Finally, 143 blood samples collected from crew members who had positive results in the first SARS-CoV-2 genetic test were analyzed in this study. A total of 269 stored blood samples were collected from patients at Nagasaki University Hospital between November 2014 and August 2019. After deduplication, 174 blood samples from patients were used as negative controls in this study (Fig 1).

The blood samples were stored at – 80°C until further antibody testing. Initially, anti-SARS-CoV-2 antibodies were detected using LFA and ECLIA by clinical technologists at Nagasaki University Hospital according to the manufacturer's instruction. Subsequently, to confirm which antibodies were influenced on LFA- and ECLIA-false-negative results in crewmember samples, IgM and IgG antibody titers against anti-SARS-CoV-2 spike-protein (SP) and nucleocapsid protein (NP) were measured using ELISA at Nagasaki University Hospital by Cellspect Co., Ltd, researchers to confirm which antibodies were influenced on LFA- and ECLIA-false-negative results. The definition of false-negative in this study is that the antibody testing was negative in the patients with negative results in the first SARS-CoV-2 genetic testing. Anti-SARS-CoV-2 antibodies were detected without notification of each result. Results from the follow-up SARS-CoV-2 genetic testing using real-time reverse transcription polymerase chain reaction RT-PCR [7] testing were obtained from the database. The follow-up genetic testing was performed on nasopharyngeal samples collected from cruise ship crew members at approximately the same time as blood samples were collected, 3 weeks after the first SARS--CoV-2 genetic testing.



Fig 1. Participant flow diagram. SARS-CoV-2, severe acute respiratory syndrome coronavirus 2; LFA, lateral flow immunochromatographic assay; ECLIA, electrochemiluminescence immunoassay; ELISA, enzyme-linked immunosorbent assay.

https://doi.org/10.1371/journal.pone.0257452.g001

Detection of anti-SARS-CoV-2 antibodies

SARS-CoV-2 IgM and IgG were detected using the SARS-CoV-2 Antibody IgM/IgG LFA testing kit (RF-NC001 and RF-NC002, Kurabo Industries, Ltd., Osaka, Japan) according to the manufacturer's instructions. The results from tests using these kits were assessed by two clinical technologists. Anti-SARS-CoV-2 total antibodies were detected by Elecsys Anti-SARS--CoV-2 (Roche Diagnostics K.K., Tokyo, Japan) using a Cobas e801 analytical unit for immunoassay tests (Roche Diagnostics K.K., Tokyo, Japan) according to the manufacturer's instructions. The judgment criteria for ECLIA were as follows: positive, cutoff index (COI) \geq 1.0; negative, COI < 1.0. To measure anti-SARS-CoV-2 IgM and IgG antibodies against SP (QuaResearch COVID-19 human IgM/IgG ELISA kit, RCOEL961N, Cellspect Co., Ltd., Iwate, Japan) and those against NP (QuaResearch COVID-19 human IgM/IgG ELISA kit, RCOEL961, Cellspect Co., Ltd., Iwate, Japan), ELISA was performed by Cellspect Co., Ltd. The antigen proteins immobilized in the ELISA kit were recombinant SP (16-1213 AA) of SARS-CoV-2 expressed in mammalian cells and recombinant NP (1-419 AA) of SARS-CoV-2 expressed in Escherichia coli. To measure IgM and IgG antibodies against SP, serum or plasma samples were diluted 1:1000 in 1% bovine serum albumin/phosphate-buffered saline with Tween detergent (PBST); to measure IgM and IgG antibodies against NP, serum or plasma samples were diluted 1:1000 in 2% non-fat milk/PBST. The absorbance was read at 450 nm using an automatic ELISA system (QRC5LB925, Cellspect Co., Ltd., Iwate, Japan) according to the manufacturer's instructions. The cutoff value determined in a previous study was used to interpret the ELISA results [8]. The cutoff values for anti-SARS-CoV-2 IgM and IgG antibody titers against SP were set at 0.25 and 0.26 and those against NP at 0.4 and 0.7, respectively [8].

Statistical analysis

All statistical analyses were performed using EZR version 1.53 (Saitama Medical Center, Jichi Medical University, Saitama, Japan), which is a graphical user interface for R (version 4.0.3; R Foundation for Statistical Computing, Vienna, Austria). Fisher's exact test or the chi-square test was used to compare categorical variables and the Mann–Whitney U test to compare continuous variables. The statistical significance level was set at P <0.05. The sensitivity (Se), specificity (Sp), positive predictive value (PPV), and negative predictive value (NPV) of the antibody testing were calculated with 95% confidence intervals (95% CI). Spearman's rank-order correlation was used to evaluate the correlation between anti-SARS-CoV-2 antibody titers in ELISA and COIs in ECLIA.

Ethics

This retrospective study was approved by the Ethics Committee of Nagasaki University Hospital (20052101) and was registered in the UMIN Clinical Trials Registry (UMIN000040402). Naso-pharyngeal samples and blood samples were anonymized and individually numbered when they were collected from the cruise ship crew members. They were sent to Nagasaki University for SARS-CoV-2 RT-PCR and antibody testing as routine care. Blood samples for the negative controls were stored anonymously at the Department of Laboratory Medicine, Nagasaki University Hospital. All data and samples were fully anonymized before we accessed them and the Ethics Committee of Nagasaki University Hospital waived the requirement for informed consent.

Results

Performance of LFA and ECLIA in detecting anti-SARS-CoV-2 antibodies

Anti-SARS-CoV-2 IgM and IgG antibodies were detected using LFA in 33 (23.1%) and 69 (48.3%) blood samples collected from crew members, respectively (Fig 2A and 2B). The blood-

(A)	(B)	(B)						
		Crew members	Neg-ctr				Crew members	Neg-ctr
1 5 4 4 -	Pos	33	0	154 (1-)	~	Pos	69	2
LFA (IGI	Neg	110	174	LFA (ig	a)	Neg	74	172
Se	e 0.231 (95%Cl, 0.164 – 0.309)				0.483 (95%Cl, 0.398 - 0.568)			
Sp	1.000 (95%)	Sp	0.989 (95%Cl, 0.959 – 0.999)					
PPV	1.000 (95%)	PPV	0.972 (95%Cl, 0.902 – 0.997)					
NPV	0.613 (95%Cl, 0.553 - 0.670)			NPV	0.601 (95%Cl, 0.638 - 0.756)			
(C)								
		Crew members	Neg-ctr					
ECLIA	Pos	112	0					
	Neg	31	174					
Se	0.783 (95%Cl, 0.707 – 0.848)							
Sp	1.000 (95%Cl, 0.969 – 1.000)							

Fig 2. Performance of LFA and ECLIA in detecting anti-SARS-CoV-2 antibodies. Anti-SARS-CoV-2 IgM (A) and IgG (B) antibodies were detected using a lateral flow immunochromatographic assay (LFA), and anti-SARS-CoV-2 total antibodies were detected using an electrochemiluminescence immunoassay (ECLIA). Blood samples were collected from cruise ship crew members who had positive results from the first SARS-CoV-2 genetic testing 3–4 weeks before blood collection. Neg-Ctr, blood samples for negative controls that were collected from patients at Nagasaki University Hospital between November 2014 and August 2019; Se, sensitivity; Sp, specificity; PPV, positive predictive value; NPV, negative predictive value.

https://doi.org/10.1371/journal.pone.0257452.g002

1.000 (95%Cl, 0.952 - 1.000)

0.848 (95%Cl, 0.792 - 0.895)

PPV

NPV

sample positivity rate for LFA-detected anti-SARS-CoV-2 IgG antibodies was significantly higher than that for IgM (p < 0.001). Anti-SARS-CoV-2 IgG antibodies were detected in all blood samples, with positive results for LFA-detected IgM antibodies. There were two positive results for anti-SARS-CoV-2 IgG antibodies in the negative controls. The Se, Sp, PPV, and NPV of LFA-detected anti-SARS-CoV-2 IgM antibodies were 0.231, 1.000, 1.000, and 0.613, respectively (Fig 2A). The Se, Sp, PPV, and NPV of LFA-detected anti-SARS-CoV-2 IgG antibodies were 0.483, 0.989, 0.972, and 0.601, respectively (Fig 2B).

Anti-SARS-CoV-2 total antibodies were detected using ECLIA in 112 (78.3%) blood samples collected from the crew (Fig 2C). There was no positive result for ECLIA-detected anti-SARS-CoV-2 total antibodies in the negative controls. Anti-SARS-CoV-2 total antibodies were detected using ECLIA in all blood samples, with positive IgG antibody results detected by LFA. The blood-sample positivity rate for ECLIA-detected anti-SARS-CoV-2 total antibodies was significantly higher than that for LFA-detected IgM and IgG antibodies (p < 0.001). The Se, Sp, PPV, and NPV of ECLIA-detected anti-SARS-CoV-2 total antibodies were 0.783, 1.000, 1.000, and 0.848, respectively (Fig 2C).

Comparison of anti-SARS-CoV-2 antibody titers between samples with false-negative and those with positive LFA results

Anti-SARS-CoV-2 SP IgM and IgG and anti-SARS-CoV-2 NP IgM and IgG antibody titers in blood samples collected from crew members were measured using ELISA (S1 Table). The median anti-SARS-CoV-2 IgM and IgG antibody titers against SP were 0.10 (0.07–0.47) and 0.27 (0.07–0.71), and those against NP were 0.18 (0.10–1.08) and 0.56 (0.18–2.82), respectively. Antibody titers were compared between blood samples with false-negative results and those with positive LFA results (Fig 3). Antibody titers for all antibodies were significantly higher in blood samples with positive results than in those with false-negative results (Fig 3A–3D). COI



Fig 3. Comparison of anti-SARS-CoV-2 antibody titers between samples with false-negative and those with positive LFA results. Anti-SARS-CoV-2 SP IgM (A) and IgG (B) and anti-SARS-CoV-2 NP IgM (C) and IgG (D) antibody titers in blood samples collected from 143 cruise ship crew members who had positive first SARS-CoV-2 genetic test results were measured using ELISA. Antibody titers were significantly lower in blood samples with false-negative LFA results (n = 74) than in those with positive LFA results (n = 69). The cutoff index in ECLIA was also significantly lower in blood samples with false-negative LFA results than in those with positive LFA results. SP, spike protein; NP, nucleocapsid protein; ELISA, enzyme-linked immunosorbent assay; LFA, lateral flow immunochromatographic assay; ECLIA, electrochemiluminescence immunoassay.

https://doi.org/10.1371/journal.pone.0257452.g003

in ECLIA was also significantly higher in blood samples with positive results than in those with false-negative results (Fig 3E). The positivity rates for anti-SARS-CoV-2 IgG antibodies against SP and NP were significantly lower in samples with false-negative results than in those with positive results (both p<0.001) (S3 Table).

Comparison of anti-SARS-CoV-2 antibody titers between samples with false-negative and those with positive ECLIA results

Antibody titers were compared between blood samples with false-negative results and those with positive ECLIA results (Fig 4). Antibody titers for all antibodies were significantly higher in blood samples with positive results than in those with false-negative results (Fig 4A–4D). The positivity rates for ELISA-detected anti-SARS-CoV-2 IgM and IgG antibodies against SP in samples with false-negative results were 0% (0/31) and 9.7% (3/31), and those against NP were 3.2% (1/31) and 0% (0/31), respectively. The positivity rates for ELISA-detected anti-SARS-CoV-2 IgG antibodies





https://doi.org/10.1371/journal.pone.0257452.g004

against SP and NP were significantly lower in samples with false-negative results than in those with positive results (p = 0.004 and p < 0.001, respectively) (<u>S4 Table</u>).

The correlation between anti-SARS-CoV-2 antibody titers in ELISA and COI in ECLIA among blood samples collected from crew members was also evaluated (Fig 5). The correlation coefficients between COI in ECLIA and anti-SARS-CoV-2 IgM and IgG antibodies against SP were 0.324 and 0.398, and those against NP were 0.332 and 0.812, respectively.

Comparison of LFA- and ECLIA-detected anti-SARS-CoV-2 antibody and result of follow-up genetic testing

Among 143 COVID-19 patients, 105 (73.4%) tested negative in the follow-up SARS-CoV-2 genetic testing using nasopharyngeal samples collected at approximately the same time as blood samples were collected. The positivity rates in crew members with negative results from the follow-up SARS-CoV-2 genetic testing were 22.9% (24/105), 47.6% (50/105), and 72.4% (76/105), respectively. The positive rates of LFA-detected anti-SARS-CoV-2 IgM and IgG antibodies and ECLIA-detected total antibodies in the crew with positive results from the follow-up SARS-CoV-2 genetic testing were 23.7% (9/38), 50.0% (19/38), and 94.7% (36/38), respectively. The positivity rate for ECLIA-detected anti-SARS-CoV-2 total antibodies was significantly higher in crew members with positive results from the follow-up SARS-CoV-2 genetic testing with positive results from the follow-up SARS-CoV-2 genetic testing with positive results from the follow-up SARS-CoV-2 genetic testing with positive results from the follow-up SARS-CoV-2 genetic testing with positive results from the follow-up SARS-CoV-2 genetic testing with positive results from the follow-up SARS-CoV-2 genetic testing than that in crew members with negative results (p = 0.003).

Discussion

This study revealed the performance of LFA and ECLIA in detecting anti-SARS-CoV-2 antibodies in asymptomatic or mild COVID-19 patients. Although all COVID-19 patients



Fig 5. Correlation between anti-SARS-CoV-2 antibody titers in ELISA and COI in ECLIA. The correlation between SARS-CoV-2 antibody titers in ELISA and cutoff indexes (COIs) in ECLIA in blood samples collected from 143 cruise ship crew members who had positive first SARS-CoV-2 genetic test results were evaluated using Spearman's rank-order correlation. SP, spike protein; NP, nucleocapsid protein; ELISA, enzyme-linked immunosorbent assay; ECLIA, electrochemiluminescence immunoassay.

https://doi.org/10.1371/journal.pone.0257452.g005

analyzed in this study were either asymptomatic or mild, we did not know their distribution. The previous study used the same ELISA kits revealed that titers of ELISA-detected anti-SARS-CoV-2 IgG antibodies against SP and NP were significantly lower in mild COVID-19 patients than those in severe and critical COVID-19 patients [8]. The titer of anti-SARS-CoV-2 IgG antibodies against SP and NP in mild COVID-19 patients in the previous study after 22 days from onset were approximately 0.5 and more than 2.0, respectively. Because they are higher than those in this study, 0.27and 0.56, respectively, it is expected that the percentage of asymptomatic COVID-19 patients was higher than that of mild COVID-19 patients in this study.

The sensitivities of LFA for both anti-SARS-CoV-2 IgM and IgG antibodies were low (0.231 and 0.483, respectively) in this study. LFA's sensitivity in this study was lower than that reported in two previous studies. One previous study reported IgM- and IgG-antibody positivity in all 24 COVID-19 patients [9]. In another study using blood samples collected from 12 COVID-19 patients, including two asymptomatic (16.7%) and seven mild patients (58.3%), the kit's sensitivity for IgM and IgG antibodies was 0.750 and 0.727, respectively [10]. ECLIA's sensitivity for anti-SARS-CoV-2 total antibodies was 0.783 in this study, which was higher than that of LFA. However, ECLIA's sensitivity in this study was lower than that reported in previous studies using the same ECLIA kit; the sensitivity in previous studies was 0.920–0.995 [11–13]. One study reported that ECLIA's sensitivity in all patients was 0.920–0.927; however, ECLIA's sensitivity in 6 (4.0%) asymptomatic patients and 37 (24.8%) mild patients was reported to be approximately 0.800 [11]. The results of this study demonstrated that ECLIA's sensitivity in detecting anti-SARS-CoV-2 antibodies was relatively low in asymptomatic or mild COVID-19 patients. Considering the results of this study, the sensitivities of LFA and

ECLIA in asymptomatic or mild COVID-19 patients were lower than the required clinical sensitivity in the guidelines [1].

In contrast to sensitivity, the specificities of LFA and ECLIA were very high in this study. These results were consistent with those of previous studies using the same kits [11–15]. Although there were no samples with false-positive results for LFA-detected anti-SARS-CoV-2 IgM antibodies and ECLIA-detected total antibodies, there were two samples with false-positive results for LFA-detected SARS-CoV-2 IgG antibodies. A previous study using the same LFA kit reported SARS-CoV-2 IgG antibody positivity in 57% (4/7) of patients with human common cold coronavirus pneumonia [9]. The results from the previous study indicated that the LFA kit used in this study had cross-reactivity with antibodies against human common cold coronavirus. Although the false-positivity rate was very low in this study, the LFA kit may not be suitable for screening COVID-19 patients or for COVID-19 serological surveillance in a population with a low prevalence of COVID-19 because specificity is especially important in such situations [1, 16].

Anti-SARS-CoV-2 IgM and IgG antibody titers against SP and NP were measured using ELISA in this study. All antibody titers were significantly lower in samples with false-negative results than in those with positive results in both LFA and ECLIA. Among the antibodies, anti-SARS-CoV-2 IgG antibodies against SP and NP were noticeably divided between blood samples with positive and those with false-negative results in both methods. Additionally, the anti-SARS-CoV-2 IgG antibody titer against NP in ELISA was strongly correlated with COI in ECLIA. The ECLIA kit used in this study detects total antibodies (IgG, IgA, and IgM) against NP [17, 18]; however, our data indicated that a low titer of anti-SARS-CoV-2 antibody against NP contributed to the false-negative results in asymptomatic and mild COVID-19 patients. A previous study reported that a combination of anti-SARS-CoV-2 antibodies against NP and SP increased the percentage of positive results [19]. In fact, the anti-SARS-CoV-2 antibody titer against SP in ELISA had a weaker correlation with COI in ECLIA in this study. However, anti-SARS-CoV-2 IgM and IgG positivity rates were also very low in COVID-19 patients with false-negative ECLIA results. Since all anti-SARS-CoV-2 antibody testing that was granted emergency use authorization by the United States Food and Drug Administration targeted SP and/or NP, it is currently difficult to improve the sensitivity of anti-SARS-CoV-2 antibody testing in asymptomatic or mild COVID-19 patients.

In this study, 73.4% of COVID-19 patients tested negative in the follow-up SARS-CoV-2 genetic testing using nasopharyngeal samples collected at approximately the same time as blood samples were collected. On the other hand, LFA and ECLIA detected anti-SARS-CoV-2 antibodies in 47.6% and 72.4% of patients with negative results in the follow-up SARS-CoV-2 genetic testing, respectively. These data support the recommendation in the guidelines, [1] although it is necessary to consider the possibility of false-negatives in anti-SARS-CoV-2 antibody testing.

This study has several limitations. First, we couldn't obtain the data regarding factors that affect antibody responses such as age, comorbid medical conditions, and immunocompromised state [1]. Since we collected the blood samples from crew members, the most patients were healthy adult. Second, only one kit each of LFA, ECLIA, and ELISA was used in this study. Thus, other kits could have possibly had higher sensitivity in asymptomatic and mild COVID-19 patients than those used in this study. However, since the ECLIA used in this study demonstrated good sensitivity compared to other ECLIA kits [11–13], it may be difficult to improve sensitivity by using others. Third, the time-course of the anti-SARS-CoV2 antibody testing was not evaluated in this study. This study focused on the time point that was reported as one of the best times for anti-SARS-CoV2 antibody testing [8, 10, 11]; however, it is possible that the optimal time for anti-SARS-CoV-2 antibody testing in asymptomatic or mild COVID-19 patients varies. Finally, the neutralizing antibodies were not evaluated. A previous study reported that the results of the ECLIA kit used in this study were correlated with SARS--CoV-2 neutralization; nevertheless, it also reported that the kit had poor negative-percent agreement for SARS-CoV-2 neutralization [20]. Therefore, further investigation is required to confirm whether COVID-19 patients with false-negative results have neutralizing antibodies.

Conclusions

The findings of this study corroborate the recommendations of the guidelines [1]. They suggest that the sensitivity of anti-SARS-CoV-2 antibody testing in asymptomatic or mild COVID-19 patients is lower than the required clinical sensitivity. In addition, it may be difficult to improve anti-SARS-CoV-2 antibody testing sensitivity at present because anti-SARS-CoV-2 IgM and IgG antibody titers against SP and NP were very low in COVID-19 patients with false-negative results. On the other hand, anti-SARS-CoV-2 antibody testing may be useful in limited contexts, such as in patients at 3–4 weeks after symptom onset but with negative SARS-CoV-2 genetic test results.

Supporting information

S1 Table. Results of LFA, ECLIA, and ELISA in blood samples collected from crew members.

(XLSX)

S2 Table. Reuslts of LFA and ECLIA in negative controls. (XLSX)

S3 Table. Comparison of positive rates for antibodies measured used ELISA between positive and false-negative result in LFA. (XLSX)

S4 Table. Comparison of positive rates for antibodies measured used ELISA between positive and false-negative result in ECLIA. (XLSX)

Acknowledgments

We would like to thank Editage (www.editage.com) for English language editing.

Author Contributions

- **Conceptualization:** Norihito Kaku, Hironori Sakai, Hiroo Hasegawa, Koichi Izumikawa, Shigeru Kohno, Katsunori Yanagihara.
- **Data curation:** Norihito Kaku, Fumitaka Nishimura, Yui Shigeishi, Rina Tachiki, Hironori Sakai, Daisuke Sasaki, Kenji Ota, Kei Sakamoto, Kosuke Kosai, Koya Ariyoshi, Hiroshi Mukae, Jiro Yasuda, Kouichi Morita.
- **Formal analysis:** Norihito Kaku, Fumitaka Nishimura, Yui Shigeishi, Rina Tachiki, Daisuke Sasaki, Kenji Ota, Hiroo Hasegawa, Koichi Izumikawa, Koya Ariyoshi, Hiroshi Mukae, Jiro Yasuda, Kouichi Morita.

Funding acquisition: Katsunori Yanagihara.

Investigation: Norihito Kaku, Fumitaka Nishimura, Yui Shigeishi, Rina Tachiki, Daisuke Sasaki, Kenji Ota, Kei Sakamoto, Kosuke Kosai.

Methodology: Norihito Kaku, Fumitaka Nishimura, Yui Shigeishi, Rina Tachiki, Hironori Sakai, Hiroo Hasegawa, Katsunori Yanagihara.

Project administration: Norihito Kaku, Shigeru Kohno, Katsunori Yanagihara.

- **Resources:** Norihito Kaku, Daisuke Sasaki, Kenji Ota, Kosuke Kosai, Koichi Izumikawa, Hiroshi Mukae, Jiro Yasuda, Kouichi Morita, Shigeru Kohno, Katsunori Yanagihara.
- **Supervision:** Norihito Kaku, Fumitaka Nishimura, Koichi Izumikawa, Koya Ariyoshi, Shigeru Kohno, Katsunori Yanagihara.

Visualization: Norihito Kaku.

- Writing original draft: Norihito Kaku, Katsunori Yanagihara.
- Writing review & editing: Fumitaka Nishimura, Yui Shigeishi, Rina Tachiki, Hironori Sakai, Daisuke Sasaki, Kenji Ota, Kei Sakamoto, Kosuke Kosai, Hiroo Hasegawa, Koichi Izumikawa, Koya Ariyoshi, Hiroshi Mukae, Jiro Yasuda, Kouichi Morita, Shigeru Kohno.

References

- Infectious Diseases Society of America Guidelines on the Diagnosis of COVID-19: Serologic Testing. [cited 11 Feb 2021]. Available: https://www.idsociety.org/practice-guideline/covid-19-guidelineserology/
- Cheng MP, Yansouni CP, Basta NE, Desjardins M, Kanjilal S, Paquette K, et al. Serodiagnostics for Severe Acute Respiratory Syndrome–Related Coronavirus 2. Ann Intern Med. 2020; 173: 450–460. https://doi.org/10.7326/M20-2854 PMID: 32496919
- Peirlinck M, Linka K, Sahli Costabal F, Bhattacharya J, Bendavid E, Ioannidis JPA, et al. Visualizing the invisible: The effect of asymptomatic transmission on the outbreak dynamics of COVID-19. Comput Methods Appl Mech Eng. 2020; 372. https://doi.org/10.1016/j.cma.2020.113410 PMID: 33518823
- Hansen CB, Jarlhelt I, Pérez-Alós L, Hummelshøj Landsy L, Loftager M, Rosbjerg A, et al. SARS-CoV-2 Antibody Responses Are Correlated to Disease Severity in COVID-19 Convalescent Individuals. J Immunol. 2021; 206: 109–117. https://doi.org/10.4049/jimmunol.2000898 PMID: 33208457
- Oran DP, Topol EJ. The Proportion of SARS-CoV-2 Infections That Are Asymptomatic. Ann Intern Med. 2021 [cited 17 Feb 2021]. https://doi.org/10.7326/M20-6976 PMID: 33481642
- Ota K, Yanagihara K, Sasaki D, Kaku N, Uno N, Sakamoto K, et al. Detection of SARS-CoV-2 using qRT-PCR in saliva obtained from asymptomatic or mild COVID-19 patients, comparative analysis with matched nasopharyngeal samples. PLoS One. 2021; 16: e0252964. https://doi.org/10.1371/journal. pone.0252964 PMID: 34111203
- National Institution of the Infectious Disease. Manual for the Detection of Pathogen 2019-nCoV Ver.2.6. 2020 [cited 12 Nov 2020]. Available: <u>https://www.niid.go.jp/niid/images/epi/corona/2019-nCoVmanual20200217-en.pdf</u>
- Imai K, Kitagawa Y, Tabata S, Kubota K, Nagura-Ikeda M, Matsuoka M, et al. Antibody response patterns in COVID-19 patients with different levels of disease severity—Japan. J Med Virol. 2021; jmv.26899. https://doi.org/10.1002/jmv.26899 PMID: 33620098
- Shibata S, Ishiguro T, Kobayashi Y, Koike M, Numano T, Shimizu Y, et al. High incidence of false-positive results of IgG antibody against SARS-CoV-2 with rapid immunochromatographic antibody test due to human common cold coronavirus infection. Respir Med Case Reports. 2020; 31: 101180. https://doi.org/10.1016/j.rmcr.2020.101180 PMID: 32766111
- Uwamino Y, Wakui M, Aoki W, Kurafuji T, Yanagita E, Morita M, et al. Evaluation of the usability of various rapid antibody tests in the diagnostic application for COVID-19. Ann Clin Biochem. 2021; 000456322098482. https://doi.org/10.1177/0004563220984827 PMID: 33334135
- Harritshøj LH, Gybel-Brask M, Afzal S, Kamstrup PR, Jørgensen CS, Thomsen MK, et al. Comparison of sixteen serological SARS-CoV-2 immunoassays in sixteen clinical laboratories. J Clin Microbiol. 2021 [cited 17 Feb 2021]. Available: http://jcm.asm.org/lookup/doi/10.1128/JCM.02596-20
- Ainsworth M, Andersson M, Auckland K, Baillie JK, Barnes E, Beer S, et al. Performance characteristics of five immunoassays for SARS-CoV-2: a head-to-head benchmark comparison. Lancet Infect Dis. 2020; 20: 1390–1400. https://doi.org/10.1016/S1473-3099(20)30634-4 PMID: 32979318
- Muench P, Jochum S, Wenderoth V, Ofenloch-Haehnle B, Hombach M, Strobl M, et al. Development and validation of the elecsys anti-SARS-CoV-2 immunoassay as a highly specific tool for determining

past exposure to SARS-CoV-2. J Clin Microbiol. 2020; 58: 1694–1714. https://doi.org/10.1128/JCM. 01694-20 PMID: 32747400

- Kittel M, Findeisen P, Muth M-C, Thiaucourt M, Gerhards C, Neumaier M, et al. Specificity testing by point prevalence as a simple assessment strategy using the Roche Elecsys® Anti-SARS-CoV-2 immunoassay. Int J Infect Dis. 2021 [cited 25 Feb 2021]. https://doi.org/10.1016/j.ijid.2021.02.024 PMID: 33578017
- Suda G, Ogawa K, Kimura M, Maehara O, Kitagataya T, Ohara M, et al. Time-dependent changes in the seroprevalence of COVID-19 in asymptomatic liver disease outpatients in an area in Japan undergoing a second wave of COVID-19. Hepatol Res. 2020; 50: 1196–1200. https://doi.org/10.1111/hepr. 13551 PMID: 32729953
- Interim Guidelines for COVID-19 Antibody Testing | CDC. [cited 25 Feb 2021]. Available: <u>https://www.cdc.gov/coronavirus/2019-ncov/lab/resources/antibody-tests-guidelines.html?deliveryName=USCDC_2067-DM29085</u>
- Tang MS, Case JB, Franks CE, Chen RE, Anderson NW, Henderson JP, et al. Association between SARS-CoV-2 Neutralizing Antibodies and Commercial Serological Assays. Clin Chem. 2020; 66: 1538– 1547. https://doi.org/10.1093/clinchem/hvaa211 PMID: 32894750
- Walker GJ, Naing Z, Ospina Stella A, Yeang M, Caguicla J, Ramachandran V, et al. SARS Coronavirus-2 Microneutralisation and Commercial Serological Assays Correlated Closely for Some but Not All Enzyme Immunoassays. Viruses. 2021; 13: 247. https://doi.org/10.3390/v13020247 PMID: 33557418
- Schnurra C, Reiners N, Biemann R, Kaiser T, Trawinski H, Jassoy C. Comparison of the diagnostic sensitivity of SARS-CoV-2 nucleoprotein and glycoprotein-based antibody tests. J Clin Virol. 2020; 129: 104544. https://doi.org/10.1016/j.jcv.2020.104544 PMID: 32663788
- Matsunaga H, Makino A, Kato Y, Murakami T, Yamaguchi Y, Kumanogoh A, et al. Radioligand Assay-Based Detection of Antibodies against SARS-CoV-2 in Hospital Workers Treating Patients with Severe COVID-19 in Japan. 2021 [cited 17 Jul 2021]. https://doi.org/10.3390/v13020347 PMID: 33672213